

To Whom It May Concern

Statement

Regarding the product Coronavirus (2019-nCoV) -Antigentest- produced by Beijing Hotgen Biotech Co., Ltd., the following statement are especially made:

The Coronavirus (2019-nCoV) -Antigentest- by our company detects N protein of SARS-CoV-2. While the mutation of the virus found in England (B.1.1.7 strain), South Africa (B.1.351 strain, **B.1.1.529 strain**), Brazil (P.1 strain) and India (B.1.617 strain, Delta variant B.1.617.2) are mainly on S protein. So our product is still effective to detect the antigen of novel coronavirus and the mutation will not affect the detection result of our product.

Beijing Hotgen Biotech Co., Ltd.

**November 28, 2021**